

SafeR Needles

PART NUMBERS

9010TMS

3ml ADULT IM/SC Retractable
Needle & Syringe Pack (23
Gauge)

9020TMS

5ml ADULT IM/SC Retractable
Needle & Syringe Pack (23
Gauge)

9030TMS

1ml PAEDIATRIC IM/SC
Retractable Needle & Syringe
Pack (25 Gauge)

Packs of 100



**One Click to
Protect,
One Use to
Prevent**

The SafeR® Retractable Safety Syringe is designed to protect healthcare workers and patients through its automatic, one-click needle retraction system. It functions like a standard syringe but ensures the needle retracts instantly after use, helping prevent needlestick injuries and reuse. Available in multiple sizes and materials, it suits a wide range of medical needs.

It incorporates key safety features such as Auto Disable (AD), Re-Use Prevention (RUP), and Sharp Injury Protection (SIP). These ensure the syringe can't be reused and that the needle is safely enclosed, minimizing injury risks. Clear calibrations and a lubricated, triple-bevelled needle also enhance dosing accuracy and comfort

The syringe works with a variety of SafeR® Sting Needles, including a blunt-fill option for drug prep. Its closed system and one-handed activation simplify use and improve safety. This design reduces the need for extra safety equipment and helps lower healthcare waste.



www.traumamedsolutions.com

sales@traumamedsolutions.com

+44(0) 1432 800962

SAFER: WHAT SETS IT APART



What sets the SafeR® Retractable Safety Syringe apart is its truly passive safety mechanism—activating automatically with no extra steps required by the user. Unlike many safety syringes that rely on manual engagement, SafeR® retracts the needle instantly upon injection completion, reducing user error and maximizing protection. Its combination of Auto Disable, Re-Use Prevention, and Sharp Injury Protection in a single, easy-to-use design ensures superior safety without compromising efficiency or comfort.

FEATURES



One-Click Needle Retraction

Automatically retracts the needle after injection, reducing the risk of needlestick injuries and contamination.



Built-In Safety Mechanisms (AD, RUP, SIP)

Includes Auto Disable, Re-Use Prevention, and Sharp Injury Protection to ensure safe, single-use performance.



Clear Dosing Calibrations

Easy-to-read markings help ensure accurate medication delivery and dosage control.



Triple-Bevelled, Lubricated Needle

Designed for smoother skin entry, improving patient comfort during injections.



1. Draw up medication in the standard way



2. Prepare and inject using aseptic technique (in line with institutional policy)



3. While the needle is in the patient, push in the plunger fully and the needle will automatically retract



4. The SafeR syringe can be thrown into the normal medical waste container after use, as the needle is safely enclosed in the cylinder.

